IL-6 WB Control

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


Materials Supplied, Storage and Stability

Components	Cat. No.	Content	Preparation	Storage at	Shelf life
QuickLine IL-6 Whole Blood Control 1: Native, human IL-6 in serum matrix (concentration 1)	MQL6BC1	1 vial à 1 mL	lyophilized; reconstitute with 1 mL aqua dest. 30 minutes before use	2 - 8 °C; Notice: After reconstitution, aliquot immediately the material and store at -20 °C!	Until the expiration date; after reconstitution 4 month at -20°C
QuickLine IL6 Whole Blood Control 2: Native, human IL-6 in serum matrix (concentration 2)	MQL6BC2	1 vial à 1 mL	lyophilized; reconstitute with 1 mL aqua. dest 30 minutes before use	2 - 8 °C; Notice: After reconstitution, aliquot immediately the material and store at -20 °C!	Until the expiration date; after reconstitution 4 month at -20°C

Application

The controls are designed only for internal quality control of Milenia® QuickLine IL-6 Whole Blood tests (MQL6B 1; test procedure see package insert). Test results can only be evaluated with Milenia® POCScan Reader.

ATTENTION:

The procedure calls for **25 µL (!)** sample volume per test (in contrast to 50 µl for whole blood samples).

Target Values

The listed target values are test method-specific values. Each laboratory should establish its own target values. Internal laboratory quality control must be adhering to current obligatory official regulations and directions.

Parameter	Cat. No.	Target (pg/ml)	Target Range (pg/ml) Preliminary
QuickLine IL6 WB Control 1	MQL6BC1	< 200	< 200
QuickLine IL6 WB Control 2	MQL6BC2	626	326 - 926

IL 6 WB Control

Interleukin-6 quality control (bi-level) for
Milenia® QuickLine IL-6 Whole Blood (MQL6B 1) –
Evaluation only with Milenia® POCScan Reader
English: Page 1-2

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Milenia Biotec GmbH