TNF- α amniotic Control

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

Materials Supplied, Storage and Stability

Components	Cat. No.	Conten	Preparation	Storage at	Shelf life
QuickLine TNF-α-amniotic Control 1: Recombinant, human TNF- α in serum matrix (concentration 1)	MQLTA C1	3 vials à 1 mL	lyophilized; reconstitute with 1 mL aqua dest. 30 minutes before use	2 - 8 °C; Notice: After reconstitution, aliquot immediately the material and store at -20 °C!	Until the expiration date; after reconstitution 4 month at -20°C
QuickLine TNF-α-amniotic Control 2: Recombinant, human TNF- α in serum matrix (concentration 2)	MQLTA C2	3 vials à 1 mL	lyophilized; reconstitute with 1 mL aqua dest. 30 minutes before use	2 - 8 °C; Notice: After reconstitution, aliquot immediately the material and store at -20 °C!	Until the expiration date; after reconstitution 4 month at -20°C

Application

The controls are designed only for internal quality control of Milenia® QuickLine TNF- α amniotic (MQLTA 1; test procedure see package insert). Test results can **only** be evaluated with Milenia® POCSan Reader.

The procedure calls for **100 μ L** sample volume per test.

Target Values

The values are test method-specific targets. Each laboratory should be establishing its own targets. Internal laboratory quality control must be adhering to current obligatory official regulations and directions.

Parameter	Cat. No.	Target (pg/ml)	Target Range (pg/ml) Preliminary
QuickLine TNF-α-amniotic Control 1	MQLTA C1	< 200	< 200
QuickLine TNF-α-amniotic Control 2	MQLTA C2	1157	434-1880

**TNF- α amniotic
Control**

TNF- α amniotic quality control (bi-level) for
Milenia[®] TNF- α amniotic (MQLTA 1) –
Evaluation only with Milenia[®] POCScan Reader
English: Page 1-2

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Milenia Biotec GmbH