TNF- α *ex vivo*

Lateral Flow Immunoassay for the quantitative measurement of human TNF- α in cell culture supernatants after LPS stimulation, interpretation by Milenia® POCScan Reader
English: Page 1-8

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

Materials Supplied, Storage and Stability

Components	Cat.-No.	Content	Preparation	Storage at	Shelf life
Test Unit: membrane coated with polyclonal (rabbit) anti-hTNF- α antibody; monoclonal (mouse) anti-hTNF- α antibody in gold conjugate	MQSEV	30 pieces	ready to use	2 - 8 °C Use immediately after opening the plastic bag!	in the closed plastic bag until the expiration date
Stimulation Reagent: Lipopolysaccharide (LPS)	MREVS	1 vial, 10 ng	lyophilized; at least 30 minutes before use, reconstitute the lyophilized LPS with 1 mL pyrogen-free distilled water. Mix thoroughly by gentle inversion before use. The reconstituted solution should be further diluted 1-in-20 with cell culture medium to a final concentration of 500 pg/mL (e.g., 1 mL reconstituted LPS solution and 19 mL Cell Culture Medium). Pipet 500 μ L of this solution into 1.5 mL reaction cups.	2-8 °C (lyophilized and reconstituted)	lyophilized until the expiration date; after reconstitution and dilution for four weeks at 2-8°C
Certificate of Analysis	CoA-QEVS	1			

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Materials Required

- Distilled, pyrogen free water
- Cell Culture Medium RPMI with Glutamine;
recommendation: Company Biochrom, Berlin/Germany, www.biochrom.de (Order code FG 1215)
- Pipets for 50, 100, 500 and 1000 μ L
- Closable (sterile not necessary) reaction cups for 1.5 mL; closable (sterile) plastic tubes for dilution
- Incubator for $37 \pm 1^\circ\text{C}$
- Centrifuge for 1,000 g

Specimen Collection and Preparation

For the whole blood stimulation heparinized blood should be taken according the rules of medical profession. It must be guaranteed that the blood collection system is pyrogen-free and no clotting takes place. The procedure calls for 50 μ L of blood per assay. The blood can be stored in the collection tubes at room temperature for 8 hours. Within this time the incubation of the sample should be started.

For determination of TNF- α in these stimulated culture supernatants, the sample should be immediately determined after stimulation. The procedure calls for 100 μ L sample volume per test.

Warnings and Precautions

All reagents of this test kit are strictly intended for *in vitro* diagnostic use only. Use by staff, who is specially informed and trained in methods which are carried out by use of *in vitro* diagnostics. Please adhere strictly to the sequence of pipetting steps provided in this protocol.

Sample material of patients (for example serum or plasma) normally used in laboratory determinations are always classified as potentially infectious. Samples of risk patients should be specially labelled and if necessary be handled in safety work benches (for example lamina flow bench).

Method and Test Principle

Milenia[®] QuickLine TNF- α *ex vivo* is the combination of standardized lipopolysaccharide (LPS) for the stimulation of human whole blood samples together with a lateral flow immunoassay for the quantitative measurement of human Tumor-Necrosis-Factor- α (TNF- α) in these cell culture supernatants.

Heparinized whole blood samples are incubated at 37°C with a standardized stimulation solution containing LPS in cell culture medium. By that the immune cells of the blood (especially the monocytes/macrophages) are stimulated to produce cytokines. By centrifugation a supernatant is collected, in which cytokines can be measured. In comparison to control experiments, the TNF- α amount produced can be indicative of activation or suppression of patient's immune system.

The culture supernatant sample is pipetted in the sample application of the test unit. TNF- α of the patient sample binds to a monoclonal anti-TNF- α antibody conjugated to gold particles. The TNF- α loaded gold particles diffuse through the membrane and overflow the test band (T). There, a polyclonal antibody specific for TNF- α is coated on the membrane; so gold particles were specifically bound and become visible as a red-blue band. Colour intensity is directly proportional to the concentration of TNF- α in the sample and band intensity increases during incubation time. The surplus of gold particles continues to diffuse over the test device. The conjugate specific antibodies printed as a line on the membrane (control band, C) capture the rest of gold conjugate and a well visible line develops during the incubation time.

Interpretation of results is only possible with Milenia[®] POCScan Reader, which allows an objective interpretation and documentation of the result.

The TNF- α concentration directly correlates with the colour intensity of the test band (T). The Milenia[®] POCScan Reader measures the colour reflection of the test band. The result is calculated according a programmed standard curve.


Assay Procedure

Attention:

- Do not interchange components of different lots.
- All components should be at room temperature (18 – 28 °C) before use.
- All components of these test kits, supplied as concentrate should be diluted to their final concentration at least 30 minutes prior to use. Mix well, but prevent foam formation.
- **A normal donor must be assayed for control!**

Further control experiments are recommended according to the respective objective.


1. Pipette **500 µL** of the ready prepared stimulation solution (LPS, 500 pg/mL) into 1.5 mL reaction cups.
2. Add **50 µL** (gently mixed) whole blood.
3. Mix preparation very carefully and incubate at 37 ± 1 °C for **4 hours**.
4. After incubation mix carefully and centrifuge **5 minutes** at 1,000 x g.
5. Pipette **100 µL** of the clear supernatant into the sample application of the test unit.
6. Incubate test unit **20 minutes** at room temperature.
7. Interpret result **immediately afterwards** with the Milenia® POCScan Reader (see page 6).

Test Performance with Milenia® POCScan Reader

Important notes:


- Prior performing the test read this instruction and the Milenia® POCScan Reader manual carefully! Follow the procedure step by step and notice the sequence!
- Before measurement, check the current lot of the used kit (see kit outerbox label!). If necessary import the data by RFID reader from the kit outerbox-label or fill them in manually according the certificate of analysis of the used kit (see Manual Milenia® POCScan Reader, Enter Calibration coefficients).

	Step	
1.		<p>Insert test unit in the Milenia® POCScan Reader.</p> <p>Pay attention to the position of test unit in the drawer!</p>
2.		<p>Close the drawer of the Milenia® POCScan Reader!</p> <p>Start measurement <u>immediately</u>.</p>
3.		<p>After scanning the result will show up in the display.</p> <p>The result must not be interpreted later!</p> <p>Only valid results are displayed.</p> <p>If result is not valid, the test must be repeated.</p> <p>For details see chapter “trouble shooting” in the manual of the Milenia® POCScan Reader.</p> <p>The manual can be downloaded from the website www.milenia-biotec.de</p>

Internal Quality Control

In any case, the control band (C) has to appear!

It is used only as functional control and can not be used for the interpretation of the test band. If the control line is not detectable after 20 minutes incubation time, the result is invalid! The analysis has to be repeated with a new test unit!

Interpretation of Results

During measurement the software controls, whether the results are valid. If the test is not valid, no test result will be shown (look "Internal Quality Control").

The normal donor used in the stimulation assay as a control of the experiment should fit the below listed TNF- α range.

It is recommended, that each laboratory checks the listed normal range, respectively establishes its own range. This TNF- α reference range should be regarded as guidelines only.

In General: Results should always be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

Expected values

In a preliminary reference study heparinized whole blood of normal donors (n = 21) was performed with the Milenia[®] QuickLine TNF- α ex vivo.

The following TNF- α values were established in the cell culture supernatants:

Range:	302 upto 1,886 pg/mL
Mean:	926 pg/mL
5% Percentile:	475 pg/mL
95 Percentile:	1,789 pg/mL

Assay Characteristics

Sample material:	Cell culture supernatant
Incubation time:	20 minutes at room temperature (18 - 28 °C)
Measurement range:	200 – 10,000 pg/mL
Lower detection limit:	200 pg/mL
Specificity:	Highly specific for the human TNF- α
Calibration:	Calibrated against 2 nd International Standard TNF- α NIBSC 88/786
High-Dose-Hook-Effect:	Not detectable for samples below 1 μ g/mL TNF- α

TNF- α *ex vivo*

Lateral Flow Immunoassay for the quantitative measurement of human TNF- α in culture supernatants after LPS stimulation, interpretation by Milenia® POCScan Reader
English: Page 1-8

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