Free IgE

Lateral Flow Immunoassay for the quantitative measurement of non complexed (Free) IgE in serum; interpretation with Milenia® POCScan Reader
English: Page 1-6

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


For further languages please visit: www.milenia-biotec.de or send an email to: info@milenia-biotec.de

Warnings and Precautions

- Read **these instructions** and the Milenia® **POCScan Reader Manual** carefully before performing the test. Follow them step by step!
- Store reagents at 2 - 8° C!
- Do not use the kit after expiration date!
- Do not use test unit if foil pouch is damaged!
- Protect test units from moisture!
- Compare batch code of your kit and the given batch code in your reader; if necessary: change batch code and parameters. (s. Milenia® **POCScan Reader Manual** –‘Enter Test Parameters’)
- Do not interchange components of packages of different batch codes!
- All samples, used units, pipettes, tips, and other contaminated material should be handled carefully and as potentially infectious.
- For professional users.

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Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage at	Shelf life
Free IgE Test Unit , membrane coated with anti-ligand and monoclonal (mouse) anti-IgE antibodies conjugated to gold particles	MQSFE	20	ready to use	2 - 8 °C <i>Use immediately after opening the foil pouch!</i>	in the <u>closed</u> foil pouch until expiration date
Sample Buffer , buffer contains ligand labelled hybrid antibody: (avian) IgY and human FcεRIα	MQFEB	1 vial à 2 mL (orange)	ready to use	2 - 8 °C	until expiration date
Certificate of Analysis	CoA-MQFE	1			

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Intended Use

- For *in vitro* diagnostic use.
- Milenia® QuickLine Free IgE is designed for quantitative measurement of non-complexed (Free) human IgE in serum for monitoring patients during immunotherapies with Anti-IgE antibodies.

Materials Required But Not Supplied

Pipet and tips: 10 µL (serum)

Specimen Collection and Preparation

Serum

Avoid haemolysis; separate serum from blood cells within two hours.

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Summary and Explanation

During the sensitization period of allergic reaction allergen-specific immunoglobulin E (IgE) are released by B-cells. These immunoglobulin E antibodies are bound by receptors on the surface of mast cells.

If a patient is again in contact with the allergenic substance (antigen), the allergen-specific IgE bound to the surface of mast cells will capture this antigen. This causes the release of mediator substances (such as histamine and inflammatory mediators) from the activated mast cell. The secretion of these substances causes the clinical symptoms of allergic reaction.

If allergies are not treated in time, the patient may develop an allergic asthma.

In recent years, modern medications for asthma therapy have been established. Their effect is based on the binding of circulating IgE, thus preventing the docking of IgE to the mast cell receptor.

These drugs must be dosed depending on the IgE levels of the patient and the body weight. Ideally, the dosage is adjusted as such that all circulating IgE molecules are bound to drug molecules. If, however, unbound IgE molecules are present, they are referred to as 'free' IgE (fIgE).

'Free' IgE should be determined to allow an individual therapy for the patient, and to monitor the effect of IgE-neutralizing drugs during treatment.

Method and Test Principle

The Milenia® QuickLine Free IgE is a lateral flow immunoassay for quantitative determination of non-complexed 'free' human immunoglobulin E (fIgE) in serum.

The sample is pipetted in the sample application port of the test unit and the sample buffer is added immediately.

The sample buffer contains a chimeric molecule, which contains binding epitopes of the human IgE receptor.

This ligand-labeled molecule binds to 'free' IgE molecules from patient's sample.

The sample and the sample buffer migrate through the membrane of the test. The formed complexes react with a monoclonal detection antibody, which is coated to gold particles. It binds to a different epitope of patient's IgE- from the sample, than the chimeric molecule.

When the complexes cross the test line (T), they are captured via the ligand molecule by the coated anti-ligand and become visible as a colored line.

The color intensity is directly proportional to the concentration of 'free' IgE in patient's sample and increases during incubation time.

The surplus of gold particles flow through the membrane to the control line (C). At the location of the C-line conjugate-specific antibodies are coated onto the membrane. These antibodies bind the conjugate with the gold particles. A clearly visible control line (C) appears.

Increasing or decreasing color's intensity of the control line cannot be used for any interpretation.

The Milenia® POCScan Reader evaluates the test unit colorimetrically after the incubation time.

The result displayed on the reader is calculated according to a stored standard curve.

Test Performance

Important notes:

- Do **not interchange** components of different lots!
- The test unit is sensitive to **humidity**. In case of damaged package do not use test unit!
- Before opening bring the test unit and sample buffer to room temperature (18 – 28 °C).
- Do **not** store at room temperature. Use the test unit immediately after opening!
- Prior performing the test read this instruction and the Milenia® POCScan Reader manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after adding the sample buffer to the patient sample in the sample port! The incubation time (300 sec) is controlled by the instrument.
- Before measurement compare **the current lot** of the used kit with the certificate of analysis (see: manual Milenia® POCScan Reader, page 12).

	Step	
1.		<p>Pipet 10 µL serum in the test unit.</p>
2.		<p>Add 2 drop sample buffer Note: Do not contaminate the vial with sample material</p>
3.		<p>Close the drawer of the Milenia® POCScan Reader! Start - SCAN - Mode, SCAN IN PROGRESS:</p> <ul style="list-style-type: none"> ■ Inkubation time (300 sec) ■ Scan ■ Result calculation
4.		<p>Result is displayed.</p> <p>If result is 'INVALID' check described items on flyer: "Notice INVALID/ A / 2009-10-16 For enquiries or error messages, please contact your local distributor or Customer Service at Milenia Biotec GmbH. E mail: cwillsch@milenia-biotec.de The Milenia® POCScan Reader manual can be downloaded from the website www.milenia-biotec.de</p>

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Interpretation of Results

The quantitative analysis of the Milenia® QuickLine Free IgE is performed by the Milenia® POCScan Reader. The calculated value is expressed in IU/mL.

During measurement the firmware controls validity of the results. If the test is not valid, no test result will be shown (see above "Functional Control").

Results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

Functional Control


In any case, the Control line (C) has to appear!

It is only used as functional control and cannot be used for the interpretation of the test line.


If the control line is not detectable after 5 minutes incubation time, the result is invalid!

The analysis has to be repeated with a new test unit!

Valid results

Invalid results

Assay Characteristics

Sample material:	10 µL serum
Incubation time:	5 minutes at room temperature (18 – 28 °C)
Sensitivity:	30 IU/mL
Dynamic range:	30 IU/mL to 3,000 IU/mL
Specificity:	Highly specific reaction with human immunoglobulin E

Limitation / Interferences

Precision: Interassay-CVs \leq 15% over entire dynamic range.

Interferences: Lipemia, bilirubin, and haemolysis have no interference influence on the results.

Note:

In this assay, free IgE is bound via the IgE receptor on the solid phase. In the event that the patient sample contains elevated levels of soluble IgE receptor (sCD23) molecules, these molecules can compete with the capture molecules for binding the free IgE. This leads to false decreased values in the assay. For this reason patients who show elevated levels of soluble IgE receptor (sCD23), e.g. patients suffering from B cell leukaemia, should be excluded from free IgE measurement with this test.

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Free IgE

Lateral Flow Immunoassay for the quantitative measurement of non complexed (Free) IgE, in serum;
interpretation with Milenia[®] POCScan Reader
English: Page 1-6

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