IL-6

**Lateral Flow Immunoassay for the quantitative measurement of human interleukin-6 (IL-6), interpretation by Milenia® POCScan Reader or PicoScan
English: Page 1-8**

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


Materials Supplied, Storage and Stability

Components	Cat.-No.	Content	Preparation	Storage at	Shelf life
Test Unit: membrane coated with monoclonal (mouse) anti-hIL-6 antibody; monoclonal (mouse) anti-hIL-6 antibody conjugated to gold particles	MQSL6	20	ready to use	2 - 8 °C Use immediately after opening the plastic bag!	until the expiration date
Chase Buffer: PBS-Buffer, pH 7.4	MQCB	1 vial à 2 mL	ready to use; Important notice: Only necessary for the interpretation with the PicoScan system!	2 - 8 °C	until the expiration date
Certificate of Analysis	CoA-MQL6	1	only necessary for the interpretation with the Milenia® POCScan Reader		

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Materials Required

Micropipet and pipet tips for 100 µL.

Specimen Collection and Preparation

For determination of interleukin 6 (IL-6) serum, plasma (EDTA, heparinized, citrated), cell culture supernatant and cerebrospinal fluid can be used (not whole blood). The blood sample should be centrifuged within 4 hours after collection of blood. If a storage longer than 8 hours at room temperature is intended, samples should be aliquoted and stored frozen at -20 °C. The procedure calls for 100 µL sample volume per test.

Warnings and Precautions

All reagents of this test kit are strictly intended for *in vitro* diagnostic use only. Use by staff, who is specially informed and trained in methods which are carried out by use of *in vitro* diagnostics. Please adhere strictly to the sequence of pipetting steps provided in this protocol.

Sample material of patients (for example serum or plasma) normally used in laboratory determinations are always classified as potentially infectious. Samples of risk patients should be specially labeled and if necessary be handled in safety work benches (for example lamina flow bench).

Method and Test Principle

The Milenia® QuickLine IL-6 test is a lateral flow immunoassay designed for the quantitative measurement of human interleukin-6 (IL-6). The sample is pipetted in the sample application. IL-6 of the patient sample binds to a first monoclonal anti-IL-6 antibody conjugated to gold particles. The IL-6-loaded gold particles diffuse through the membrane and overflow the test band (T). There is a second monoclonal antibody specific for IL-6 coated on the membrane; so gold particles were specifically bound and become visible as a red-blue band. Color intensity is directly proportional to the concentration of IL-6 in the sample and intensifies during incubation time. The surplus of gold particles continues to diffuse over the test device. The conjugate specific antibodies printed as a line on the membrane (control band, C) captures the rest of gold conjugate and a well visible line develops during the incubation time.



Interpretation of results is only possible with Milenia® POCScan Reader or PicoScan systems, which allow an objective interpretation and documentation of the result.

The IL-6 concentration directly correlates with the color intensity of the test band (T). The Milenia® POCScan Reader measures the color reflection of the test band, the PicoScan the color intensity of the test band respectively. The result is calculated according a programmed standard curve.

1) Test Performance with Milenia® POCScan Reader

Important notes:

- The test is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do not store at room temperature!
- Prior performing the test read this instruction and the Milenia® POCScan Reader manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after pipetting the sample into the sample port! The incubation time of 20 minutes is controlled by the instrument.
- Before measurement scan the current lot of the used kit by the RFID reader from the outerbox label of the kit or check the lot manually with the certificate of analysis of the kit (look manual Milenia® POCScan Reader).


	Step	
1.		<p>Pipet 100 µL sample in the round sample port of the test unit (see arrow in the figure).</p> <p>Attention: Pipetting of chase buffer is <u>not</u> necessary!</p>
2.		<p>Close the drawer of the Milenia® POCScan Reader!</p> <p>Start measurement <u>immediately</u>.</p>
3.		<p>After 20 minutes the result shows up in the display.</p> <p>The result must not be interpreted later!</p> <p>Only valid results are displayed.</p> <p>If result is not valid, the test must be repeated.</p> <p>For details see chapter “trouble shooting” in the manual of the Milenia® POCScan Reader.</p> <p>The manual can be downloaded from the website www.milenia-biotec.de</p>

2) Test Performance with PicoScan

Important notes:


- The test is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit and chase buffer to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do not store at room temperature!
- Prior performing the test read this instruction and the PicoScan manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after adding the chase buffer to the patient sample in the sample port! The last incubation time of 5 minutes is controlled by the instrument.

1. Pipet **100 µL** sample into the sample application of the test unit (see figure).
2. After **15 minutes** incubation **add 2 drops chase buffer** directly with the bottle in the sample application of the test unit.
3. After further **5 minutes** of incubation the test result is interpreted using the **PicoScan system**.

The result is automatically saved and can be directly read from the computer screen.

The results should not be interpreted afterwards!

Interpretation of Results

During measurement the software controls, whether the results are valid. If the test is not valid, no test result will be shown (look "Internal Quality Control").

You can find further informations in the PicoScan manual, which you can download from our website www.milenia-biotec.de.



In General: Positive results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

Internal Quality Control



In any case, the control band (C) has to appear!

It is used only as functional control and can not be used for the interpretation of the test band. If the control line is not detectable after 20 minutes incubation time, the result is invalid! The analysis has to be repeated with a new test unit!

Valid results

Invalid results

Assay Characteristics

Sample material:	Serum, plasma (EDTA, heparinized, citrated), culture supernatant, cerebrospinal fluid
Incubation time:	20 minutes at room temperature (18 - 28 °C)
Specificity:	Highly specific for human IL-6
Measurement range:	50 – 10,000 pg/mL
Sensitivity:	50 pg/mL
Reference values:	< 100 pg/mL in healthy blood donors
High-Dose-Hook-Effect:	Not detectable for samples \leq 500,000 pg/mL IL-6
Interference Factors:	Lipemia / Bilirubin / Hemolysis no influence for values > 100 pg/mL if necessary, control values \leq 100 pg/mL visually
Calibration:	Calibrated with NIBSC WHO standard 89/548 (PicoScan, NIBSC IL-6)
Method Comparison:	60 plasma samples of ICU patients and 20 serum samples from healthy donors were analyzed; Milenia [®] QuickLine IL-6 result were independently interpreted by two operators; 77.7 % (62/80) show identical results in comparison to IMMULITE [®] IL-6 (DPC, Los Angeles, USA)

IL-6

Lateral Flow Immunoassay for the quantitative measurement of human interleukin-6 (IL-6), interpretation by Milenia® POCScan Reader or PicoScan
English: Page 1-8

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