## IL-6 Whole Blood

**Lateral Flow Immunoassay for the quantitative measurement of human interleukin-6 (IL-6) in whole blood, interpretation by Milenia® POCScan Reader or PicoScan**  
English: Page 1-8

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


## Materials Supplied, Storage and Stability

Components	Cat.-No.	Content	Preparation	Storage at	Shelf life
<b>Test Unit:</b> Membrane coated with an anti-ligand; monoclonal (mouse) anti-hIL-6 antibody in gold conjugate	MQSL6B	20	ready to use	2 - 8 °C Use <b>immediately</b> after opening the plastic bag!	in the closed plastic bag until the expiration date
<b>IL-6 Whole Blood Sample Buffer, lot specific:</b> Buffered solution with a ligand coupled to a monoclonal (mouse) anti-hIL-6 antibody	MQSB	1 vial à 1.5 ml (pink cap)	ready to use	2 - 8 °C	until the expiration date
<b>Transfer pipettes 50 µL, disposables</b>	MQP50	20	see test procedure		
<b>Certificate of Analysis</b>	CoA-MQL6B	1	only necessary for the interpretation with the Milenia® POCScan Reader		

Material Safety Data Sheets are available on request (look as well [www.milenia-biotec.de](http://www.milenia-biotec.de)).

## Specimen Collection and Preparation

For determination of interleukin 6 (IL-6) capillary, venous and cord blood can be used. Make sure that the sample is not clotted (preferably blood donation with anticoagulants like heparin, EDTA; no citrate). The blood sample should be determined immediately after blood donation. The procedure calls for 50 µL whole blood sample volume per test.

## Warnings and Precautions

All reagents of this test kit are strictly intended for *in vitro* diagnostic use only. Use by staff, who is specially informed and trained in methods which are carried out by use of *in vitro* diagnostics. Please adhere strictly to the sequence of pipetting steps provided in this protocol.

Sample material of patients (for example serum or plasma) normally used in laboratory determinations are always classified as potentially infectious. Samples of risk patients should be specially labeled and if necessary be handled in safety work benches (for example lamina flow bench).

## Method and Test Principle

The Milenia<sup>®</sup> QuickLine IL-6 Whole Blood test is a lateral flow immunoassay designed for the quantitative measurement of human interleukin-6 (IL-6) in whole blood. The sample is pipetted in the sample application. Sample buffer containing a first monoclonal anti-IL-6 antibody coupled to a ligand is added. IL-6 of the patient sample binds to this antibody and to a second monoclonal anti-IL-6 antibody conjugated to gold particles in the application pad. Erythrocytes and all other cells are prevented from further flowing into the evaluation window by a special membrane pad. The IL-6 loaded gold particles diffuse through the membrane and overflow the test band (T). There is an anti-ligand coated on the membrane; so IL-6 loaded gold particles were specifically bound and become visible as a red-blue band. Color intensity is directly proportional to the concentration of IL-6 in the sample and intensifies during incubation time. The surplus of gold particles continues to diffuse over the test device. The conjugate specific antibodies printed as a line on the membrane (control band, C) captures the rest of gold conjugate and a well visible line develops during the incubation time.



Interpretation of results is only possible with Milenia<sup>®</sup> POCScan Reader or PicoScan systems, which allow an objective interpretation and documentation of the result.

The IL-6 concentration directly correlates with the color intensity of the test band (T). The Milenia<sup>®</sup> POCScan Reader measures the color reflection of the test band, the PicoScan the color intensity of the test band respectively. The result is calculated according a programmed standard curve.

## 1) Test Performance with Milenia® POCScan Reader

### Important notes:

- Do **not interchange** components of different lots!
- The test is sensitive to **humidity**. In case of damaged package do not use test unit!
- Before opening bring the test unit and sample buffer to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do **not** store at room temperature!
- Prior performing the test read this instruction and the Milenia® POCScan Reader manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after adding the sample buffer to the patient sample in the sample port! The incubation time of 25 minutes is controlled by the instrument.
- Before measurement scan **the current lot** of the used kit by the RFID reader from the outerbox label of the kit or check the lot manually with the certificate of analysis of the kit (look manual Milenia® POCScan Reader).

	Step	
1.		<p><b>Pipet 50 µL patient blood in the test unit.</b></p> <p>Handle the provided transfer pipet like the following: For taking up the blood sample hold the pipette on the blood drop <u>without pushing</u>. The blood is automatically filled in upto the black mark (see upper arrow in the figure).</p> <p>Put the transfer pipet with the patient sample directly into the round sample port of the test unit (see lower arrow in the figure).</p> <p>Press fingers on the pinhole, located at the black label of the transfer pipet, squeeze the pipet bulb until the blood has been completely released into the sample port.</p>
2.		<p><b>Add <u>immediately afterwards</u> 1 drop sample buffer</b> directly from the vial to the blood in the same round sample port of the test unit (see arrow in the figure).</p>
3.		<p>Close the drawer of the Milenia® POCScan Reader!</p> <p><b>Start measurement <u>immediately</u>.</b></p>
4.		<p>After <b>25 minutes</b> the <b>result</b> shows up in the display.</p> <p><b>The result must not be interpreted later!</b></p> <p><b>Only valid results are displayed.</b></p> <p>If result is not valid, the test must be repeated.</p> <p>For details see chapter “trouble shooting” in the manual of the Milenia® POCScan Reader.</p> <p>The manual can be downloaded from the website <a href="http://www.milenia-biotec.de">www.milenia-biotec.de</a></p>

## 2) Test Performance with PicoScan

### Important notes:

- **Do not interchange** components of different lots!
- The test is sensitive to **humidity**. In case of damaged package do not use test unit!
- Before opening bring the test unit and sample buffer to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do **not** store at room temperature!
- Prior performing the test read this instruction and the PicoScan manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after adding the sample buffer to the patient sample in the sample port!

1. Apply **50 µL** whole blood on the sample application bottom of the test unit (see figures).

The provided transfer pipette is used in the following way:

- For taking up the blood sample hold the pipette on the blood drop **without pushing**. The blood is automatically filled in upto the black mark.
- For giving off the blood sample hold closed the small hole at the black mark, put the pipette directly on the sample application and squeeze the bulb.

2. **Directly afterwards** add **only 1 drop** sample buffer directly with the bottle in the sample application of the test unit.

3. Subsequently put at once the test unit into the PicoScan System and start the automatic measurement („Interleukin-6 WB 25 min“).

The test unit will be automatically measured after an incubation time of **25 minutes** in the system.

The result is automatically saved and can be directly read from the computer screen.

### **The results should not be interpreted afterwards!**

Software interface showing patient information and test results:

Name: SMITH, ANNA    DOB: 16.08.2007    ID:    Next QC measure

Patient Sample: ?    Operator: Administrator

Test	Date	Time
IL-6 WB	17.08.2007	13:14:27
QC1	17.08.2007	12:49:05

Test: Interleukin-6  
Conc. pg/ml: 1358  
Range: 1000-2000 pg/ml

Interleukin\_6\_WB\_25min v5 UK Test valid

Comment:

## Interpretation of Results

During measurement the software controls, whether the results are valid. If the test is not valid, no test result will be shown (look "Internal Quality Control").

You can find further informations in the PicoScan manual, which you can download from our website [www.milenia-biotec.de](http://www.milenia-biotec.de).


**In General:** Positive results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

## Internal Quality Control


**In any case, the control band (C) has to appear!**

It is used only as functional control and can not be used for the interpretation of the test band. If the control line is not detectable after 25 minutes incubation time, the result is invalid! The analysis has to be repeated with a new test unit!

Valid results

Invalid results

## Assay Characteristics

Sample material:	Whole blood like capillary blood, cord blood, venous blood (heparin, EDTA, <u>no citrate</u> )
Incubation time:	25 minutes at room temperature (18 - 28 °C)
Measurement range:	50 – 10,000 pg/mL
Sensitivity:	50 pg/mL
Specificity:	highly specific for human IL-6
Reference values:	< 100 pg/mL in healthy blood donors
Calibration:	calibrated against Milenia® QuickLine IL-6 (MQL6 1)
High-Dose-Hook-Effect:	not detectable for samples $\leq 500,000$ pg/mL IL-6
Interference Factors:	Lipemia, bilirubin and hemolysis have no interference influence on the results

## IL-6 Whole Blood

Lateral Flow Immunoassay for the quantitative measurement of human interleukin-6 (IL-6) in whole blood, interpretation by Milenia® POCScan Reader or PicoScan  
English: Page 1-8

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The screenshot shows the PicoScan Reader 490 Secure software interface. The patient information includes Name: SMITH, ANNA, DOB: 16.08.2007, and ID. The test results table shows:

Test	Date	Time
IL-6 WB	17.08.2007	13:14:27
QC1	17.08.2007	12:49:05

The visual representation of the test strip shows two lanes: Control and IL-6. The Control lane shows a single line, and the IL-6 lane shows two lines. The test result is Interleukin-6, with a concentration of 1358 pg/ml and a range of 1000-2000 pg/ml. The test is valid.

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