TNF- α amniotic

Lateral Flow Immunoassay for the quantitative measurement of human TNF- α in amniotic fluid, interpretation by Milenia[®] POCScan Reader
English: Page 1-8

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

Materials Supplied, Storage and Stability

Components	Cat.-No.	Content	Preparation	Storage at	Shelf life
Test Unit: membrane coated with polyclonal (rabbit) anti-hTNF- α antibody; monoclonal (mouse) anti-hTNF- α antibody in gold conjugate	MQSTA	20	ready to use	2 - 8 °C Use immediately after opening the plastic bag!	in the closed plastic bag until the expiration date
Transfer pipettes 100 μL, disposables	MQP100	20	see test procedure		
Certificate of Analysis	CoA-MQLTA	1			

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Materials Required

Optionally: Micropipet and pipet tips for 100 μ L; the disposable transfer pipet provided with the kit is identical in use.

Specimen Collection and Preparation

For determination of TNF- α only amniotic fluid can be used. The sample should be immediately determined after donation. The procedure calls for 100 μ L sample volume per test.

Warnings and Precautions

All reagents of this test kit are strictly intended for *in vitro* diagnostic use only. Use by staff, who is specially informed and trained in methods which are carried out by use of *in vitro* diagnostics. Please adhere strictly to the sequence of pipetting steps provided in this protocol.

Sample material of patients (for example serum or plasma) normally used in laboratory determinations are always classified as potentially infectious. Samples of risk patients should be specially labeled and if necessary be handled in safety work benches (for example lamina flow bench).

Method and Test Principle

The Milenia[®] QuickLine TNF- α amniotic test is a lateral flow immunoassay designed for the quantitative measurement of human TNF- α in amniotic fluid. The sample is pipetted in the sample application. TNF- α of the patient sample binds to a monoclonal anti-TNF- α antibody conjugated to gold particles. The TNF- α loaded gold particles diffuse through the membrane and overflow the test band (T). There, a polyclonal antibody specific for TNF- α is coated on the membrane; so gold particles were specifically bound and become visible as a red-blue band. Colour intensity is directly proportional to the concentration of TNF- α in the sample and band intensity increases during incubation time. The surplus of gold particles continues to diffuse over the test device. The conjugate specific antibodies printed as a line on the membrane (control band, C) capture the rest of gold conjugate and a well visible line develops during the incubation time



Interpretation of results is only possible with Milenia[®] POCScan Reader, which allow an objective interpretation and documentation of the result.

The TNF- α concentration directly correlates with the color intensity of the test band (T). The Milenia[®] POCScan Reader measures the color reflection of the test band. The result is calculated according a programmed standard curve.

Test Performance with Milenia® POCScan Reader

Important notes:

- The test is sensitive to humidity. In case of damaged package do not use test unit!
- Before opening bring the test unit to room temperature (18 – 28 °C). Use the test unit immediately after opening. Do not store at room temperature!
- Prior performing the test read this instruction and the Milenia® POCScan Reader manual carefully! Follow the procedure step by step and notice the sequence!
- Start measurement immediately after pipetting the sample into the sample port! The incubation time of 20 minutes is controlled by the instrument.
- Before measurement scan the current lot of the used kit by the RFID reader from the outerbox label of the kit or check the lot manually with the certificate of analysis of the kit (look manual Milenia® POCScan Reader).


	Step	
1.		<p>Pipet 100 µL amniotic fluid in the test unit.</p> <p>The provided (capillary) transfer pipette is used in the following way:</p> <ul style="list-style-type: none"> - Close the black marked pinhole (see upper arrow) and use it like a normal squeeze pipette. - Take the sample up only to the black mark! - Put the pipette directly into the sample application of the test unit (see lower arrow), squeeze the bulb and dispense in this way the complete fluid volume in the sample application of the test unit.
2.		<p>Close the drawer of the Milenia® POCScan Reader!</p> <p>Start measurement <u>immediately</u>.</p>
3.		<p>After 20 minutes the result shows up in the display.</p> <p>The result must not be interpreted later!</p> <p>Only valid results are displayed.</p> <p>If result is not valid, the test must be repeated.</p> <p>For details see chapter “trouble shooting” in the manual of the Milenia® POCScan Reader.</p> <p>The manual can be downloaded from the website www.milenia-biotec.de</p>

Internal Quality Control


In any case, the control band (C) has to appear!

It is used only as functional control and can not be used for the interpretation of the test band. If the control line is not detectable after 20 minutes incubation time, the result is invalid! The analysis has to be repeated with a new test unit!

Valid results

Invalid results

Interpretation of Results

During measurement the software controls, whether the results are valid. If the test is not valid, no test result will be shown (look "Internal Quality Control").

In General: Positive results should be verified concerning the entire clinical status of the patient. Additionally every decision for therapy should be taken individually.

Assay Characteristics

Sample material:	Amniotic fluid
Incubation time:	20 minutes at room temperature (18 - 28 °C)
Measurement range:	200 – 10,000 pg/mL
Lower detection limit:	200 pg/mL
Specificity:	Highly specific for the human TNF- α
Calibration:	Calibrated against 2 nd International Standard TNF- α NIBSC 88/786
High-Dose-Hook-Effect:	Not detectable for samples below 1 μ g/mL TNF- α
Interference Factors:	Lipemia and bilirubin have no interference influence on the results High haemolysis can cause false negative results; therefore control visually, if necessary

TNF- α amniotic

Lateral Flow Immunoassay for the quantitative measurement of human TNF- α in amniotic fluid, interpretation by Milenia® POCScan Reader
English: Page 1-8

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