## Total IgE

**Lateral Flow Immunoassay for the quantitative measurement of Total IgE in combination with Milenia® POCScan Reader**  
English: Page 1-12

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

## Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage at	Shelf life
<b>Total IgE Test Cassette</b> , membrane coated with anti-ligand and monoclonal (mouse) anti-IgE antibodies conjugated to gold particles	MQSTE	20	ready to use	2 - 8 °C Use <b>immediately</b> after opening the foil bag!	in the sealed foil bag until expiration date
<b>Assay Buffer</b> , buffer containing ligand labeled polyclonal (goat) anti-IgE antibodies	MQTEB	1 vial à 1.5 mL	ready to use	2 - 8 °C	until expiration date
<b>Lancets</b> , (sterile)	3090	20	see Test Performance	room temperature	until expiration date
<b>Transfer pipettes</b> , (calibrated for 20 µL)	MQP20	1 package à 20 pipettes	see Test Performance	room temperature	until expiration date
<b>Cleansing swabs</b>	3104	20	ready to use	room temperature	until expiration date
<b>Adhesive tapes</b>	3105	20	ready to use	room temperature	until expiration date
<b>Manual for test and reader</b>	MQTE	1			
<b>Certificate of Analysis</b>	CoA-MQTE	1			

Material Safety Data Sheets are available on request (look as well [www.milenia-biotec.de](http://www.milenia-biotec.de)).

## Specimen Collection and Preparation

The determination of Total IgE calls for 20 µL capillary blood.

## Warnings and Precautions

All reagents of this test kit are strictly intended for *in vitro* diagnostic use only. Use by staff, who is specially informed and trained in methods which are carried out by use of *in vitro* diagnostics. Please adhere strictly to the steps of the work flow provided in this protocol.

Sample material of patients (for example blood) normally used in laboratory determinations are always classified as potentially infectious. Samples of risk patients should be specially labeled and if necessary be handled in safety work benches (for example lamina flow bench).

## Method and Test Principle


In combination with the Milenia<sup>®</sup> POCScan Reader the Milenia<sup>®</sup> QuickLine Total IgE is a lateral flow immunoassay which allows the quantitative determination of total immunoglobulin E (IgE) from whole blood.

Immediately after pipetting the whole blood sample to sample port of the test cassette, one drop of assay buffer is added. IgE antibodies of the patient sample then specifically bind to the ligand-labeled polyclonal anti-IgE antibodies present in the assay buffer. These immune complexes migrate into the test strip and there associate with the gold conjugate, which are coated with monoclonal anti-IgE antibodies. The gold particles overflow the test line (T) where an anti-ligand is immobilized. Only the immune complexes consisting of ligand labeled anti-IgE antibodies, patient's IgE and the gold conjugate will bind to the anti-ligand. Over time a reddish test band is generated. Gold conjugate-specific antibodies, which are printed as a line on the membrane (control line, C), catch unbound gold conjugate; so a clearly visible line is generated during the test time.


Interpretation of results is only possible with the Milenia<sup>®</sup> POCScan Reader, which allows an objective interpretation and documentation of the result.

The total IgE concentration directly correlates with the color intensity of the test band (T). The Milenia<sup>®</sup> POCScan Reader measures the color reflection of the test band and calculates the result according a programmed standard curve.

## Milenia® POCScan Reader Overview

### Battery position and replacement:













- Please use Nickel-Metal-Hybride (NiMH) batteries only, which offer longer durability.
- **ATTENTION: Recharging of batteries is working only, if POCScan reader is switched ON!**




## Test Performance



### Important notes:




- The test is sensitive to humidity. **In the case of damaged package do not use test unit!**
- Before opening bring test cassette and assay buffer to room temperature (18 – 28 °C). Use the test cassette immediately after opening. Do not use kit reagents from different lots!
- **Please read the instruction carefully prior performing the test; follow the procedure step by step and notice the sequence !**
- **Start measurement immediately after adding the assay buffer to the patient sample in the sample port!** The incubation time of 7 minutes is controlled by the instrument.
- For identification the operating elements of Milenia® POCScan Reader see overview on page 5.

<b>A) PREPARATIONS</b>		
	<b>Step</b>	* Remark: All figures (fig.) are example displays only!
<b>1.</b>		<p>Turn Milenia® POCScan Reader on. For that, press “ON/OFF” button  for two/three seconds.</p> <p>A self test/ self calibration runs automatically. Menu „SCAN“ is automatically selected.</p>
<b>2.</b>		<p>Take the test cassette out of the aluminium foil bag.</p> <p>Open the drawer of Milenia® POCScan Reader and place the test cassette into the holder.</p>
<b>3.</b>		<p>Start menu “SCAN”.</p> <p>For that, press “ENTER” button .</p>
<b>4.</b>		<p>Enter USER ID and SAMPLE ID with respective buttons ● (UP/DOWN/BACK/FORWARD/ENTER).</p> <p>Select “Next” and press .</p>

	Step	Continuation (Preparations)
5.		<p>Compare lot number and product code displayed by the Milenia® POCScan Reader, with the lot number and the product code of the outer box label of the kit.</p> <p>If they are identical, select “Yes” and press </p> <p>The display automatically switches to <b>step 10</b>.</p>
6.		<p>If the displayed lot number is <u>not</u> identical with the lot number of the used kit, lot number and respective curve parameters need to be changed.</p> <p>In this case select “No” and press </p> <p>The Menu “Enter coefficients” opens.</p>
7.		<p>Select Lot, Prm A, Prm F and Prm G with respective buttons  and change all respective data according the values listed in the “<b>Certificate of Analysis</b>” of the currently used Milenia® QuickLine Total IgE kit.</p> <p>Select “Next” and press </p>
8.		<p>Select Prm H, Prm J, and Prm K with respective buttons  and change all respective data according the values listed in the “<b>Certificate of Analysis</b>” of the currently used Milenia® QuickLine Total IgE kit.</p> <p>Select “Next” and press  The display switches back to the Scan Mode</p>
9.		<p>Compare displayed Checksum with the Checksum listed in the “<b>Certificate of Analysis</b>” of the currently used Milenia® QuickLine Total IgE kit.</p> <p>If the checksums are identical, select “Yes” and press </p>
9a.		<p>If Checksum is <u>not</u> identical, please select “No”, press  and accurately check in <b>step 7 and 8</b> the different parameters of the “<b>Certificate of Analysis</b>”.</p>
10.		<p>The menu “START” appears.</p> <p><b>Do <u>not</u> yet activate this menu!</b></p> <p><b>Now draw blood from the finger tip of the patient (look part B, BLOOD DONATION).</b></p>

<b>B) BLOOD DONATION</b>		
	<b>Step</b>	
1.		Take the lancet and push the cap into the body until it clicks. Hold the lancet as described in the figure.
2.		Twist off and discard the lancet cap. Do not touch the release button (see white square in the fig.). The lancet is now ready to use.
3.		The marked areas of patient's fingers are ideal for the prick. Clean the area of the selected finger with the cleansing swab and do not touch the finger until you prick. Wait until the finger is dry.
4.		Press the platform of the lancet firmly against the chosen site on the finger and press the release button.
5.		Massage the finger from the hand to the direction of the prick (see arrow in the fig.), until a drop of blood is formed.
6.		Hold the transfer pipet from the side to the blood drop. Massage the finger until the transfer pipet is automatically filled up to the black label of the pipet (see arrow in the fig.).  <b>Pipet immediately the blood sample into the test cassette (look part C, TEST).</b>
7.		If needed, fix an adhesive tape over the prick of patient's finger.

<b>C) TEST</b>	
<b>Step</b>	
<b>1.</b>	 <p><b>Pipet the patient blood into the test cassette.</b></p> <p>For that, handle the transfer pipet like the following:</p> <p>Put the transfer pipet with the patient sample directly into the round sample port of the test cassette (see lower arrow in the fig.).</p> <p>Press fingers on the pinhole, located at the black label of the transfer pipet (see upper arrow in the fig.), squeeze the pipet bulb until the blood has been completely released into the sample port.</p>
<b>2.</b>	 <p><b>Add <u>immediately afterwards</u> 1 drop assay buffer</b> directly from the vial to the blood in the same round sample port of the test cassette (see arrow in the fig.).</p>
<b>3.</b>	 <p><b>Close the drawer of the Milenia® POCScan Reader!</b></p> <p><b>Start measurement (look part D, MEASUREMENT).</b></p>

D) MEASUREMENT		
	Step	
1.		<p>Now start <b>immediately</b> measurement by pressing "Enter" button .</p>
2.		<p>Incubation time of 420 seconds (7 min) is counting down (test measurement runs then automatically).</p>
3.		<p><b>After 7 minutes the result of Milenia® QuickLine Total IgE shows up in the display. The result must not be interpreted later!</b></p> <p>"Next test" is selected automatically; press  to return to "Main Menu".</p> <p>Now, the system is ready for a new measurement, <u>or</u></p> <p>switch off the system by pressing the "ON/OFF" button  for two/three seconds.</p>

## Interpretation of Results

The interpretation of Milenia® QuickLine Total IgE can be done only in combination with the Milenia® POCScan Reader. The calculated value is expressed in U/mL.

## Internal Quality Control

Only valid results are shown in the display of Milenia® POCScan Reader.

If the result is shown up as invalide, the analysis has to be repeated with a new test cassette!

## Assay Characteristics

Sample material:	20 µL capillary blood
Incubation time:	7 minutes at room temperature (18 – 28 °C)
Sensitivity:	15 U/mL
Dynamic range:	30 U/mL to 2,000 U/mL
Specificity:	Highly specific reaction with human immunoglobulin E
Method Comparison:	High correlation with a world wide established laboratory method for determination of Total IgE from serum/plasma
Precision:	Interassay-CVs ≤ 15% over entire dynamic range
Interferences:	Lipemia, bilirubin and hemolysis have no interference influence on the results

**Total IgE**

Lateral Flow Immunoassay for the quantitative measurement of Total IgE in combination with Milenia® POCScan Reader  
English: Page 1-12

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