IL-6 Control

**Interleukin-6 quality control (bi-level) for
Milenia® QuickLine IL-6 (MQL6 1)
– Evaluation only with Milenia® POCScan Reader
English: Page 1-2**

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


Materials Supplied, Storage and Stability

Components	Cat. No.	Content	Preparation	Storage at	Shelf life
QuickLine IL-6 Control 1: Native, human IL-6 in serum matrix (concentration 1)	MQL6C1	3 vials à 1 mL	lyophilized; reconstitute with 1 mL aqua dest. 30 minutes before use	2 - 8 °C; Notice: After reconstitution, aliquot immediately the material and store at -20 °C!	Until expiration date; after reconstitution 4 month at -20°C
QuickLine IL-6 Control 2: Native, human IL-6 in serum matrix (concentration 2)	MQL6C2	3 vials à 1 mL	lyophilized; reconstitute with 1 mL aqua dest. 30 minutes before use	2 - 8 °C; Notice: After reconstitution, aliquot immediately the material and store at -20 °C!	Until expiration date; after reconstitution 4 month at -20°C

Application

The controls are designed only for internal quality control of Milenia® QuickLine IL-6 tests (MQL6 1; test procedure see package insert). Test results can **only** be evaluated with Milenia® POCScan Reader.

The procedure calls for **100 µL** sample volume per test.

Target Values

The values are test method-specific targets. Each laboratory should be establishing its own targets. Internal laboratory quality control must be adhering to current obligatory official regulations and directions.

Parameter	Cat. No.	Target (pg/ml)	2s-Range (pg/ml)	3s-Range (pg/ml)
QuickLine IL-6 Control 1	MQL6C1	< 50	< 50	< 50
QuickLine IL-6 Control 2	MQL6C2	467	259-675	155-779

IL-6 Control

Interleukin-6 quality control (bi-level) for
 Milenia® QuickLine IL-6 (MQL6 1)
 – Evaluation only with Milenia® POCScan Reader
 English: Page 1-2

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Milenia Biotec GmbH