



Milenia® QuickLine Total IgE



Informations about Milenia Biotec:

- ▶ Located in Gießen/Germany
- ▶ Founded in 2000
- ▶ Development of rapid tests for the diagnosis of inflammatory and allergic diseases
- ▶ Manufacturer of CE labeled products according to IVD Directive 98/79 EC
- ▶ Certified according to EN ISO 13485:2003+AC:2007, EN ISO 9001:2000

▶ **Milenia Biotec GmbH**
Versailler Str. 1
D-35394 Gießen
Germany

Phone: ++49 641 94888 3-0
Fax: ++49 641 94888 3-80

info@milenia-biotec.de
www.milenia-biotec.de

Order Informations:

Parameter	Order Code	Pack Size
▶ Milenia® QuickLine Total IgE	MQTE 1	20 tests
▶ Milenia® POCScan Reader	MSCAN 1	1 system

- ▶ *If you are interested in these products, please fill in this page and fax it to: **++49 641 948883 - 80***
- ▶ Yes, I am interested in testing your products, please contact me

▼ *Title, Surname, First name*

Institution

Department

Street

City, Postcode

Tel. / Fax

E-mail

Date, Sign



Quantitative Total IgE results from whole blood within 7 minutes

Milenia® QuickLine Total IgE

- ▶ Rapid test for quantitative determination of total IgE from whole blood
- ▶ Kit box contains: lancets, transfer pipets, alcohol swabs, adhesive tapes, test cassettes, assay buffer
- ▶ Material supplied for 20 determinations
- ▶ Calibrated with 2nd WHO standard
- ▶ High correlation with a widely used laboratory method
- ▶ Package insert in your mother tongue
- ▶ RFID-coded method on outer box label
- ▶ CE-marked IVD

Product code: MQTE 1



Features of the Milenia® QuickLine/POCScan-System

- ▶ Rapid test performance (7 minutes)
- ▶ Measurement of total IgE directly from whole blood
- ▶ Dynamic range: 30 – 2,000 U/mL
- ▶ Only valid results are reported
- ▶ Internal QC and validity check points
- ▶ RFID coding of all lot-specific features
- ▶ Easy handling
- ▶ Automatic report of the total IgE level:

Example Result

Scan Mode
tIgE: 247 U/ml Test ID: 0001
<input type="button" value="NEXT TEST"/> <input type="button" value="DEL"/> <input type="button" value="PRINT"/>



- ▼ Prick finger tip



- ▼ Pick up calibrated whole blood volume

- ▼ Transfer blood sample to test cassette

- ▼ Add one drop of the assay buffer



- ▼ Insert test cassette into Milenia® POCScan Reader

- ▼ Start measurement

- ▼ Read test result after 7 minutes

Milenia® POCScan Reader

- ▶ Colorimetric scanner for reading rapid tests
- ▶ Stand alone usage; no PC required
- ▶ Dimensions (in cm): 18 x 17,5 x 4,5
- ▶ Weight: 200 g
- ▶ External RFID reader
- ▶ USB interface
- ▶ External power supply
- ▶ Battery powered
- ▶ Internal battery charger
- ▶ Protective plastic case
- ▶ CE-marked instrument

Product code: MSCAN 1